Instruction

WIESLAB[®] ASCA IgG semi quantitative kit

Enzyme immunoassay for detection of IgG autoantibodies
against the mannan structure of *Saccharomyces cerevisiae*

Break apart microtitration strips (12x8) 96 wells
Store the kit at +2-8° C
For in vitro diagnostic use only


Document No. LABEL-DOC-0039 v2.0
December, 2018

WIESLAB[®] ASCA IgG semi quantitative kit

English:	page	2
Français:	page	12
Español:	página	19
Deutsch:	Seite	26
Italiano:	pagina	33
Dansk:	side	43
Norsk:	side	50
Svenska:	sida	57

ASCA 151

INTENDED USE

The Wieslab® ASCA test kit is an enzyme-linked immunosorbent assay (ELISA) for detection and semi-quantitation of IgG antibodies to mannan of *Saccharomyces cerevisiae* in human sera. The assay is used to detect antibodies in a single serum specimen. The results of the assay are to be used as an aid to the diagnosis of Crohn's disease.
FOR *IN VITRO* DIAGNOSTIC USE.

Summary and explanation

Antibodies to the mannan structure of *Saccharomyces cerevisiae* has been known for several years (1) and their association to Crohn's disease discussed (2, 3). Their diagnostic use have, however, not been appreciated until the ANCA test was established and it was shown that ANCA occurred in ulcerative colitis. Thus, it became possible to use these two tests in combination to differentiate Crohn's disease from UC (4,5). The clinical presentation of these patients are often similar and then diagnosis depend on clinical, radiographic, endoscopic and laboratory data and especially in pediatric patients a serological test can be of use to avoid biopsy. The exact sensitivity and specificity is not yet settled. Studies indicate that 40-80% of UC patients have ANCA and 0-20% of Crohn's patients have ANCA, while 0-10% of UC patients have ASCA and 50-70% of Crohn's patients have ASCA. In one study the combination of positive ANCA with negative ASCA for diagnosis of UC had a sensitivity of 57% and specificity of 97%, while the combination of negative ANCA with positive ASCA for the diagnosis of Crohn's had a sensitivity of 49% and a specificity of 97% (6).

Principle of the Wieslab® ASCA assay

The wells of the microtitre strips are coated with mannan structure of *Saccharomyces cerevisiae*. During the first incubation, specific antibodies in diluted serum, will bind to the antigen coating.

The wells are then washed to remove unbound antibodies and other components.

A conjugate of alkaline phosphatase-labelled antibodies to human IgG binds to the antibodies in the wells in this second incubation.

After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. The amount of bound antibodies correlates to the colour intensity and is measured in terms of absorbance (optical density (OD)). The absorbance is then calculated against a calibrator curve and the results are given in arbitrary U/mL.

Warnings and precautions

- For in vitro diagnostic use.
- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The concentrations of ASCA in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity.

- Material safety data sheets for all hazardous components contained in this kit are available on request from Svar Life Science.

BUF	WASH	30X
DIL		
CONJ		
CAL		

CONTROL	+
CONTROL	-
SUBS	pNPP

Warning

Contains ProClin 300:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

- H317: May cause an allergic skin reaction.
 P264: Wash hands thoroughly after handling.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P302+352: IF ON SKIN: Wash with plenty of soap and water.
 P333+313: If skin irritation or rash occurs: Get medical advice/attention.

Specimen collection

The ASCA assay is for use with serum. Handle as if capable of transmitting infectious agents.

Avoid using sera which are icteric, lipemic and hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used. Store serum between 2-8° C if testing will take place within five days.

If specimens are to be kept for longer periods, store at -20° C or colder. Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results. The NCCLS provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

Kit components and storage of reagents

- One frame with break apart wells (12x8) coated with mannan structure of *Saccharomyces cerevisiae* and one lid sealed in a foil pack with a dry pack.
 - 1.5 mL negative control (NC) containing human serum in diluent.
 - 1.5 mL positive control (PC) containing human serum in diluent.
 - 13 mL conjugate containing alkaline phosphatase-labelled antibodies to human IgG (colour blue).
 - 2 x 32 mL Diluent (Dil) containing PBS (red colour).
 - 13 mL Substrate pNPP-solution.
 - 30 mL wash solution 30x concentrated.
 - Five calibrators containing human serum in diluent. 1.5 mL Cal 1 = 160 U/mL, 1.5 mL Cal 2 = 80 U/mL, 1.5 mL Cal 3 = 40 U/mL, 1.5 mL Cal 4 = 20 U/mL, 1.5 mL Cal 5 = 10 U/mL.
- All reagents in the kit, except the wash solution, are ready for use and should be stored at 2-8° C. Remove only the number of strips needed for testing, resealing the aluminium package carefully.

Materials or equipment required but not provided

Microplate reader with filter 405 nm.

Precision pipettes with disposable tips.

Washer for strips, absorbent tissue, tubes and a timer.

PROCEDURE

All solutions should be used at room temperature.

Remove only the number of wells needed for testing, resealing the aluminium package carefully.

Incubate in all steps at room temperature (20-25° C) with lid. Incubate serum for 30 minutes, conjugate for 30 minutes and substrate for 60 minutes (± 10 minutes).

Preparation of washing solution

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37°C water bath until the crystals have dissolved before dilution of wash solution. Dilute 10 mL of the 30x concentrated wash solution in 290 mL distilled water. When stored at 2-8° C, the diluted wash solution is stable until the date of expiration of the kit.

Dilution of serum and incubation

Dilute the patient's serum 1/100 with diluent (495 µL diluent +5 µL serum).

Quantitative assay: Pipet 100 µL/well in duplicate of diluent (as a blank), Calibrator 1, 2, 3, 4, 5, NC, PC and diluted patient's serum (P) according to the diagram below.

Qualitative assay: Pipet 100 µL/well in duplicate Calibrator 3, NC, PC and diluted patient's serum (P) according to the diagram below. Incubate for 30 minutes.

Quantitative assay:

	1	2	3	4	5
A	Dil	Cal 4	P1		
B	Dil	Cal 4	P1		
C	Cal 1	Cal 5	P2		
D	Cal 1	Cal 5	P2		
E	Cal 2	NC	etc		
F	Cal 2	NC			
G	Cal 3	PC			
H	Cal 3	PC			

Qualitative assay:

1	2	3	4	5
Cal 3	P2			
Cal 3	P2			
NC	etc.			
NC				
PC				
PC				
P1				
P1				

After serum incubation

Wash 3 times with 300 µL washing solution / well, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

Adding conjugate

Add 100 µL conjugate to each well. Incubate for 30 minutes.

After conjugate incubation

Wash as before.

Adding substrate solution

Add 100 µL substrate solution to each well, incubate for 60 minutes (± 10 minutes).

Read the absorbance at 405 nm on a microplate reader.

Calculations

Quantitative assay:


Subtract the OD value for the blank from the other OD values.

Construct a calibrator curve by plotting the OD against the unit values of the 5 calibrators.

The five calibrators provided have been assigned arbitrary values of 160 U/mL for calibrator 1, 80 U/mL for calibrator 2, 40 U/mL for calibrator 3, 20 U/mL for calibrator 4 and 10 U/mL for calibrator 5. Read the unit value of the patient from the constructed curve. Values greater than 160 should be reported as >160, or reassay them with a higher dilution.

Arbitrary U/mL have been adopted, as no generally recognised international standard exists for expressing ASCA titres.

Example:	Calibrator	U/mL	Absorbance
	1	160	2.11
	2	80	1.78
	3	40	1.41
	4	20	0.97
	5	10	0.58

A sample with an absorbance value of 1,37 will read on the X-axis as having 39 U/mL of ASCA

Important: The curve is an example and should not be used for actual patient interpretation.

Qualitative assay:

An optical density (OD) ratio for PC, NC each patient sample is calculated as follows:

$$\text{OD ratio} = \frac{\text{OD patient sample}}{\text{OD of Cal 3}}$$

The patient sample is negative if the OD ratio is < 0.8 and positive if the OD ratio is >0.9.

Quality Control

Quantitative assay:

The OD for the negative control should be less than that of calibrator 4.

The OD for calibrator 1 should be > 1.4

The value for the positive control, see lot certificate.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cut-off. It is recommended to assay an additional control at the assay cut-off.

Qualitative assay:

The OD for calibrator 3 should be between 1.0-1.8.

The ratio for the negative control should be < 0.8

The ratio for the positive control see lot certificate.

If any of the values are not within their respective ranges, the test should be considered invalid and the test should be repeated. Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organisations. Refer to NCCLS C24-A for guidance on appropriate QC practices.

Interpretation of results

Quantitative assay	Qualitative assay
< 29 U/mL = Negative	OD ratio < 0.8 = Negative
29-35 U/mL = Equivocal ; Retest, if still equivocal retest by an alternative method or test a new sample	0.8-0.9 = Equivocal
>35 U/mL = Positive	> 0.9 = Positive

An equivocal or positive sample should always be confirmed by a quantitative assay.

Limitations

The individual patient's antibody titre can not be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardisation of results.

The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used in combination with clinical symptoms and the results of other available tests.

Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive for ASCA antibodies with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies. Therapy should not be started on basis of a positive ASCA result. Initiation or changes in treatment should not be based on changes in ASCA concentration alone, but rather on careful clinical observation.

Expected results

A few percent of normal healthy individuals are found to be positive in ASCA. The exact sensitivity and specificity is not yet settled. Studies indicate that 0-10% of UC patients have ASCA and 50-70% of Crohn's patients have ASCA.

Performance characteristics

Table 1. Clinical sensitivity and specificity. A total of 417 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total number	Negative <29U/mL	Equivocal 29-35U/mL	Positive >35U/mL
Blood donors:	72	72	0	0
WG, MP, SLE, RA:	40	37	2	1
Celiac:	48	45	0	3
Ulcerative Colitis:	140	137	0	3
Crohn's disease:	117	75	3	39

WG = Wegener's granulomatosis, MP = microscopic polyangiitis, SLE = systemic lupus erythematosus, RA = rheumatoid arthritis

Clinical sensitivity (Equivocal samples excluded from calculations)

Crohn's disease = 39/114 = 34 % 95% CI = 25.5-42.9 %

Clinical specificity (Equivocal samples excluded from calculations)

WG, MP, RA, SLE	= 37/38 = 97 %	95% CI = 86.2-99.9 %
Celiac	= 45/48 = 94 %	95% CI = 82.8-98.7 %
Ulcerative Colitis	= 137/140 = 98 %	95% CI = 93.9-99.6 %
Normal	= 72/72 = 100 %	95% CI = 95.0-100 %

The 95% confidence interval (CI) was calculated using the exact method.

Table 2. Inter-assay precision was determined by testing three different samples in duplicate. Results were obtained for four different runs.

Sample	Mean value U/mL	SD	CV %
PK	40	3	8
K157	3	6	
K355	8	14	

Table 3. Intra-assay precision was determined by testing one sample in 22 wells.

Sample	Mean value U/mL	SD	CV %
1	39	2	5

Table 4. Linearity. The values were determined for serial two-fold dilutions of three positive sera. The values were compared to log 2 of dilution by standard linear regression. The data indicates that the assay has a linear relationship with serum dilution.

Serum	neat	1:2	1:4	1:8	1:16	r
1	162	52	29	13	4	0.985
2	84	42	24	12	3	0.998
3	116	47	24	11	3	0.996
4	107	48	24	12	3	0.999
5	188	69	36	18	7	0.991

Troubleshooting

Problem:	Possible causes:	Solution:
Control values out of range.	<ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Improper dilution. 4. Optical pathway not clean. 	<ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Repeat test. 4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.
All test results negative.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Antigen coated plate inactive. 	<ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagents. Repeat test. 2. Check for obvious moisture in unused wells. Wipe bottom and reread.
All test results yellow.	<ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated. 3. Improper dilution of serum. 	<ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water solution used to prepare. 3. Repeat test.
Poor precision.	<ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5%. 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells. 6. Improper pipetting. 	<ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips.

References:

1. **Main J, McKenzie H, Kerr MA, Robson D, Pennington CR, Parrat D.** Antibody to *Saccharomyces cerevisiae* (bakers yeast) in Crohns disease. *BMJ* 1988, 297, 1105-1106.
2. **Barnes R, Allan S, Taylor-Robinson C, Finn R, Johnson P.** Serum antibodies reactive with *Saccharomyces cerevisiae* in inflammatory bowel disease: Is IgA antibody a marker for Crohns disease? *Int Arch Allergy Appl Immunol* 1990, 92, 9-15.
3. **Sendid B, Quinton J, Charrier G, Goulet O, Cortot A, Grandbastien B, Poulain D, Colobel J.** Anti-*Saccharomyces cerevisiae* mannan antibodies in familial Crohns disease. *Am J Gastroenterol* 1998, 93, 1306-10.
4. **Quinton J, Sendid B, Reumax D, Duthilleul P, Cortot A, Grandbastien B, Charrier G, Targan S, Colombel J, Poulain D.** Anti-*Saccharomyces cerevisiae* mannan antibodies combined with antineutrophil cytoplasmic autoantibodies in inflammatory bowel disease: prevalence and diagnostic role. *GUT* 1998, 42, 788-91.
5. **Kugathasan S, Werlin S.** Measurement of p-ANCA and ASCA in screening for IBD in young children. *Inflamm Bowel Dis* 1999, 5, 283-84.
6. **Vermeire S, Peeters M, Rutgeerts P.** Diagnostic approach to IBD. *Hepatogastroenterology* 2000, 47, 44-48.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Forklaring til symboler. Symbolforklaring. Förklaringar till symboler.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partnummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuil de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
 96	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Fortynning. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Konjugat. Conjugato. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Vaskebuffert 30 x konc. Vaskeløsning 30x kons. Tvättbuffert 30x konc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
CAL X	Calibrator. Etalon. Calibrador. Kalibrator. Calibratore. Kalibrator. Kalibrator.
CONTROL X	Control. Contrôle. Control. Kontrolle. Controllo. Kontrol. Kontroll. Kontroll.

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References:

1. **Main J, McKenzie H, Kerr MA, Robson D, Pennington CR, Parrat D.** Antibody to *Saccharomyces cerevisiae* (bakers yeast) in Crohns disease. *BMJ* 1988, 297, 1105-1106.
2. **Barnes R, Allan S, Taylor-Robinson C, Finn R, Johnson P.** Serum antibodies reactive with *Saccharomyces cerevisiae* in inflammatory bowel disease: Is IgA antibody a marker for Crohns disease? *Int Arch Allergy Appl Immunol* 1990, 92, 9-15.
3. **Sendid B, Quinton J, Charrier G, Goulet O, Cortot A, Grandbastien B, Poulain D, Colobel J.** Anti-*Saccharomyces cerevisiae* mannan antibodies in familial Crohns disease. *Am J Gastroenterol* 1998, 93, 1306-10.
4. **Quinton J, Sendid B, Reumax D, Duthilleul P, Cortot A, Grandbastien B, Charrier G, Targan S, Colombel J, Poulain D.** Anti- *Saccharomyces cerevisiae* mannan antibodies combined with antineutrophil cytoplasmic autoantibodies in inflammatory bowel disease: prevalence and diagnostic role. *GUT* 1998, 42, 788-91.
5. **Kugathasan S, Werlin S.** Measurement of p-ANCA and ASCA in screening for IBD in young children. *Inflamm Bowel Dis* 1999, 5, 283-84.
6. **Vermeire S, Peeters M, Rutgeerts P.** Diagnostic approach to IBD. *Hepatogastroenterology* 2000, 47, 44-48.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Forklaring til symboler. Symbolforklaring. Förklaringar till symboler.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partnummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuil de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Fortynning. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Konjugat. Conjugato. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Vaskebuffert 30 x konc. Vaskeløsning 30x kons. Tvättbuffert 30x konc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
CAL X	Calibrator. Etalon. Calibrador. Kalibrator. Calibratore. Kalibrator. Kalibrator.
CONTROL X	Control. Contrôle. Control. Kontrolle. Controllo. Kontrol. Kontroll. Kontroll.

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Advarsel

Indeholder ProClin 300:

Blanding af: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

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SVAR LIFE SCIENCE AB