Instruction

WIESLAB® Capture MPO-ANCA

Enzyme immunoassay for detection of autoantibodies
against myeloperoxidase

Microtitration 96 wells
Store the kit at +2-8° C
For in vitro diagnostic use only

Document No. LABEL-DOC-0048, 2.0
December, 2018

WIESLAB® Capture MPO-ANCA

English:	page	2
Français:	page	13
Español:	page	20
Deutsch:	seite.....	27
Italiano:	pagina	38
Português:	página	45
Dansk:	side.....	52
Norsk:	side.....	59
Svenska:	sida.....	66

INTENDED USE

The Wieslab® Capture MPO-ANCA test kit is an enzyme-linked immunosorbent assay (ELISA) for detection and quantitation of IgG antibodies to myeloperoxidase (MPO) in human sera. The results of the assay are to be used as an aid to the diagnosis of microscopic polyangiitis (MP). The assay is intended for use in patients with signs and symptoms consistent with MP. It is not intended for screening a healthy population. The analysis should be performed by trained laboratory professionals. FOR IN VITRO DIAGNOSTIC USE.

Summary and explanation

ANCAs (anti-neutrophil cytoplasmic antibodies) are a family of autoantibodies related to vasculitis and inflammatory disorders. When c-ANCA was shown to relate to Wegener's granulomatosis (WG), interest in ANCAs has increased steadily, and today these antibodies are considered to be major diagnostic tools for the investigation of systemic vasculitis.

The first method to detect ANCA was indirect immunofluorescence (IIF) performed on ethanol fixed granulocytes (1). This method yields two patterns, a cytoplasmic staining of the granulocyte denoting the presence of c-ANCAs, and a perinuclear staining denoting the presence of p-ANCAs. IIF was followed by ELISAs using purified proteins (2, 3).

The granulocyte is full of granules each with many different proteins. It was shown that antibodies from systemic vasculitis patients bind to the alpha fraction containing the azurophil granules. The most important proteins were proteinase 3 (PR3) (4, 5) and myeloperoxidase (MPO)(6). PR3 is a serine protease with a molecular weight of 29kD, and MPO is a dimer with a molecular weight of 140kD. Thus antibodies to proteinase 3 are termed PR3-ANCA, and antibodies to myeloperoxidase are termed MPO-ANCA.

Approximately 80-90% of WG patients manifest PR3-ANCA and 5-15% MPO-ANCA. One category of vasculitis is microscopic polyangiitis (MP). Most patients with active MP are characterised by positive ANCA test results, MPO-ANCA being more frequent than PR3-ANCA (7-12).

An international workgroup has developed an international standard for MPO-ANCA serology. The Wieslab® Capture MPO-ANCA is standardized against the AF-CDC international standard (Reference Human Serum 15, code IS2720).

Principle of the Wieslab® Capture MPO-ANCA ELISA

The wells of the microtitre plate are coated with purified anti-MPO monoclonal antibody and myeloperoxidase. During the first incubation, specific antibodies in diluted serum, will bind to the antigen coating.

The wells are then washed to remove unbound antibodies and other components.

A conjugate of alkaline phosphatase-labelled antibodies to human IgG binds to the antibodies in the wells in this second incubation.

After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. The amount of bound antibodies correlates to the colour intensity and is measured in terms of absorbance (optical density (OD)). The absorbance is then calculated against a calibrator curve and the results are given in IU/mL adapted to the AF-CDC international standard for MPO.

Warnings and precautions

- For in vitro diagnostic use.
- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommend that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient samples to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The concentrations of anti-MPO in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity.
- Material safety data sheets for all hazardous components contained in this kit are available on request from Svar Life Science.

BUF	WASH	30X
DIL		
CONJ		
CAL		

CONTROL	+
CONTROL	-
SOLN	STOP
SUBS	pNPP

Warning

Contains ProClin 300:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

- H317: May cause an allergic skin reaction.
 P264: Wash hands thoroughly after handling.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P302+352: IF ON SKIN: Wash with plenty of soap and water.
 P333+313: If skin irritation or rash occurs: Get medical advice/attention.

Specimen collection

The Wieslab® Capture MPO-ANCA assay is for use with serum. Handle as if capable of transmitting infectious agents.

Avoid using sera, which are icteric, lipemic and hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used.

Store serum between 2-8° C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20° C or colder. Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.

The CLSI provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

Kit components and storage of reagents

- One frame with 96 wells coated with monoclonal anti-myeloperoxidase/myeloperoxidase with lid, sealed in a foil pack with a dry pack.
- 1.5 mL negative control (NC) containing human serum in diluent.
- 1.5 mL positive control (PC) containing human serum in diluent.
- 13 mL conjugate containing alkaline phosphatase-labelled antibodies to human IgG (blue colour).
- 2 x 32 mL Diluent (Dil) containing PBS (red colour).
- 13 mL Substrate pNPP.
- 13 mL Stop Solution.
- 30 mL wash solution 30x concentrated.
- Six calibrators (five calibrators, Cal 2-6, containing human serum) in diluent. 1.5 mL Cal 1 = 0 IU/mL, 1.5 mL Cal 2 = 2 IU/mL, 1.5 mL Cal 3 = 10 IU/mL, 1.5 mL Cal 4 = 30 IU/mL, 1.5 mL Cal 5 = 100 IU/mL, 1.5 mL Cal 6 = 200 IU/mL.

All reagents in the kit are ready for use except wash solution and should be stored at 2-8° C. Remove only the number of wells needed for testing, resealing the aluminium package carefully.

Materials or equipment required but not provided

- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- Washer for strips, absorbent tissue, tubes and a timer.

PROCEDURE

All solutions should be used at room temperature. Do not open the foil packaged plate until it has reached room temperature as damp from condensation may have a negative effect on the antigen. Incubate in all steps at room temperature (18-25° C) with lid. Incubate serum for 60 minutes, conjugate for 30 minutes and substrate for 30 minutes.

Preparation of washing solution

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37°C water bath until the crystals have dissolved before dilution of wash solution. Dilute 10 mL of the 30x concentrated wash solution in 290 mL distilled water. When stored at 2-8° C, the diluted wash solution is stable until the date of expiration of the kit.

Dilution of serum and incubation

Dilute the patient's serum 1/100 with diluent (990 µL diluent +10 µL serum).

Pipet 100µL/well in duplicate of Calibrator 1, 2, 3, 4, 5, 6, PC, NC and diluted patient's serum (P) according to the diagram below. Incubate for 60 minutes.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Cal 1	Cal 5	P1									
B	Cal 1	Cal 5	P1									
C	Cal 2	Cal 6	P2									
D	Cal 2	Cal 6	P2									
E	Cal 3	PC	etc									
F	Cal 3	PC										
G	Cal 4	NC										
H	Cal 4	NC										

After serum incubation

Wash 3 times with 300 µL washing solution / well, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

Adding conjugate

Add 100 µL conjugate to each well. Incubate for 30 minutes.

After conjugate incubation

Wash as before.

Adding substrate solution

Add 100 µL substrate pNPP to each well, incubate for 30 minutes.

Adding stop solution


Add 100 µL stop solution to each well and read the absorbance at 405 nm on a microtiterplate reader within 2 hours.

Calculations

Construct a calibrator curve by plotting the OD against the IU/mL values of the 6 calibrators. The six calibrators provided have been titrated to match the AF-CDC standard resulting in the values of 0 IU/mL for calibrator 1, 2 IU/mL for calibrator 2, 10 IU/mL for calibrator 3, 30 IU/mL for calibrator 4, 100 IU/mL for calibrator 5 and 200 IU/mL for calibrator 6. Read the IU/mL value of the patient from the constructed curve. Values greater than 200 should be reported as >200, or reassay them with a higher dilution. The Wieslab® Capture MPO-ANCA is standardized against the AF-CDC international standard.

Example:	Calibrator	IU/mL	Absorbance
	1	0	0.100
	2	2	0.224
	3	10	0.606
	4	30	1.093
	5	100	1.555
	6	200	1.845

A sample with an absorbance value of 1.437 will read on the X-axis as having 65 IU/mL of MPO-ANCA. In this example a 4 parameter logistic curve fit has been applied.

Important: The curve is an example and should not be used for actual patient interpretation.

Quality Control

The OD for calibrator 1 should be < 0.2.

The OD for calibrator 6 should be > 1.0.

The value for the positive and negative controls, see lot certificate.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cut-off. It is recommended to assay an additional control at the assay cut-off. As the positive control is ready-to-use it does not indicate an eventual dilution error by the user. It is recommended to use an internal control for this purpose.

If any of the control values are not within their respective ranges, the test should be considered invalid and the test should be repeated. Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organisations. Refer to CLSI C24-A for guidance on appropriate QC practices.

Interpretation of results

< 5 IU/mL = **Negative**

5-7 IU/mL = **Equivocal**; Retest, if still equivocal retest by an alternative method or test a new sample

>7 IU/mL = **Positive**

Limitations

The individual patient's antibody titre cannot be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardisation of results. The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used in combination with clinical symptoms and the results of other available tests. Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive for MPO antibodies with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies. Individuals receiving mouse anti-human antibodies for treatment or diagnosis, or those patients who have been otherwise exposed to mouse immunoglobulin, may produce Human Anti-Mouse Antibodies (HAMA). These antibodies can interfere with assays using mouse monoclonal antibodies and may cause falsely elevated levels. Immunosuppressive therapy should not be started on basis of a positive ANCA result. Initiation or changes in treatment should not be based on changes in ANCA concentration alone, but rather on careful clinical observation.

Expected results

MPO-ANCA are rarely found in normal healthy individuals. The Wieslab® Capture MPO-ANCA was tested with 120 normal sera. All were found to be negative. One new patient with MPO-ANCA is expected per 100000 individuals per year. (7) MPO-ANCA activity is more frequently found in microscopic polyangiitis (MP) than in Wegener's granulomatosis. However, some true vasculitis patients test negative in both IIF and ELISA. The Wieslab® Capture MPO-ANCA was tested with 54 WG sera of which 6 were found to be positive. The Wieslab® Capture MPO-ANCA was tested with 79 MP sera of which 46 were found to be positive.

Performance characteristics

Table 1. Clinical sensitivity and specificity. A total of 348 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total number	Negative <5 IU/mL	Equivocal 5-7 IU/mL	Positive >7 IU/mL
Blood donors:	120	120	0	0
WG:	54	46	2	6
MP:	79	32	1	46
GBM:	55	42	4	9
SLE:	14	14	0	0
RA:	14	13	0	1
Sjögren's syndrome	12	12	0	0

WG = Wegener's granulomatosis, MP = microscopic polyangiitis
 SLE = systemic lupus erythematosus RA = rheumatoid arthritis
 GBM = glomerular basement membrane

Clinical sensitivity (equivocal samples excluded from calculations)

WG = $6/52 = 11.5\%$ 95% CI = 0.44 – 23.4%
 MP = $46/78 = 59.0\%$ 95% CI = 47.3 – 70.0%

Clinical specificity (equivocal samples excluded from calculations)

GBM = $42/51 = 82.4\%$ 95% CI = 69.1 – 91.6%
 SLE = $14/14 = 100\%$ 95% CI = 76.8 - 100%
 RA = $13/14 = 92.9\%$ 95% CI = 66.1 – 99.8 %
 Sjögren's syndrome = $12/12 = 100\%$ 95% CI = 73.5 - 100%
 Donors = $120/120 = 100\%$ 95% CI = 97.0 - 100%

The 95% confidence interval (CI) was calculated using the exact method.

Table 2. Agreement between the Capture MPO-ANCA kit and an alternative capture MPO-ELISA.

A total of 348 frozen retrospective sera were assayed. The following table summarises the results.

Wieslab® Capture-MPO-ANCA KIT					
Alternative ELISA		Positive	Equivocal	Negative	Total
	Positive	58	5	2	65
	Equivocal	0	0	1	1
	Negative	4	2	276	282
	Total	62	7	279	348

Sera falling in the equivocal range were excluded from the following calculations

Positive Percent Agreement:	$58/60 = 96.7\%$	95% CI = 88.5– 99.6%
Negative Percent Agreement:	$276/280 = 98.6\%$	95% CI = 89.5 – 99.6%
Overall Percent Agreement:	$334/340 = 98.2\%$	95% CI = 91.4 – 99.4%

The 95% confidence interval (CI) was calculated using the exact method.

Table 3. Batch to batch variation was determined by testing six different samples in eight replicates on three different batches.

	1	2	3	4	5	6
Mean value (IU/mL)	8	6	9	33	44	58
SD	0.85	0.74	0.62	1.69	2.38	4.35
% CV	11	13	7	5	5	8

Table 4. Inter-assay precision was determined by testing six different samples in eight replicates at three separate occasions.

	1	2	3	4	5	6
Mean value (IU/mL)	8	6	10	31	44	58
SD	0.57	0.74	1.61	5.20	6.31	6.86
% CV	7	11	17	17	14	12

Table 5. Intra-assay precision was determined by testing six different samples in eight replicates at one occasion.

	1	2	3	4	5	6
Mean value (IU/mL)	7	8	11	28	36	52
SD	0.49	0.93	0.69	1.46	3.07	2.21
% CV	7	11	6	5	9	4

Table 6. Dilution recovery was determined by testing five serial dilutions for three different samples.

Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution Corrected % Recovery
1	1/100	187	187	100
	1/200	101	94	108
	1/400	54	47	116
	1/800	28	23	120
	1/1600	14	12	120
Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution Corrected % Recovery
2	1/100	87	87	100
	1/200	45	44	103
	1/400	23	22	106
	1/800	12	11	110
	1/1600	6	5	110
Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution Corrected % Recovery
3	1/100	81	81	100
	1/200	43	41	106
	1/400	21	20	104
	1/800	10	10	99
	1/1600	5	5	99

Since patient sera contain heterogeneous antibody populations some samples may exhibit non-linearity, especially at very high sample dilutions.

Troubleshooting

Problem:	Possible causes:	Corrective action:
One or more calibrators out of range.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Improper pipetting. 	<ol style="list-style-type: none"> 1. Test invalid. Repeat test. 2. Test invalid. Repeat test.
Control values out of range.	<ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Improper dilution. 4. Optical pathway not clean. 	<ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Repeat test. 4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.
All test results negative.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Antigen coated plate inactive. 	<ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagent. Repeat test. 2. Check for obvious moisture in unused wells. Wipe bottom and reread.
All test results yellow.	<ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated. 3. Improper dilution of serum. 	<ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water solution used to prepare. 3. Repeat test.
Poor precision.	<ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5%. 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells. 6. Improper pipetting. 	<ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips.

References:

1. **Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): Manual of biological markers of disease. Dordrecht, Kluwer Academic Publishers; 1993, A9:1-14.
2. **Rasmussen N, Sjolín C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). *J Immunol Methods* 1990,127:139-145.
3. **Wieslander J.** How are anti neutrophil cytoplasmic autoantibodies detected? *Am J Kidney Dis* 1991,18:154-158.
4. **Jenne DE, Tschopp J, Ludemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. *Nature* 1990, 346:520.
5. **Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
6. **Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
7. **Jennette C, Falk R.** Small vessels vasculitis. *N Engl J Med* 1997, 337:1512-1525.
8. **Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. *Kidney Int* 1995, 48: 844-850.
9. **Cohen Tervaert JW, Huitema MG, Hene RJ, Sluiter WJ, The T, van der Hem GK et al.** Prevention of relapses in Wegener's granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre. *Lancet* 1990; 336:709-711.
10. **Savigne J, Davis D, Falk DJ, Jennette JC, Wiik A.** ANCA and associated diseases; A review of clinical and laboratory features. *Kidney Int* 2000, 57:846-862.
11. **Savigne J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. *Am J Clin. Pathol* 1999, 111:507-513.
12. **Csernok E, Ahlquist D, Ullrich S, Gross W.L.** A critical evaluation of commercial immunoassays for antineutrophil cytoplasmic antibodies directed against proteinase 3 and myeloperoxidase in Wegener's granulomatosis and microscopic polyangiitis. *Rheumatol* 2002, 41:1313-1317.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklarung der Symbole. La spiegazione di simboli. Explicao dos smbolos Forklaring til symboler. Symboler som brukes p etiketter. Frklaringar till symboler.

	Batch code. Numero de lot. Nmero de lote. Chargen-Nummer. Numero di lotto. Cdigo do lote. Partinummer. Lot nummer. Satsnummer.
	Catalogue number. Rfrence catalogue. Nmero de catlogo. Katalog-Nummer. Numero di catalogo. Nmero catalogo. Katalognummer.
	Use-by date. Date de premption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udlbsdato. Utlpsdato. Anvnd fre.
	Temperature limit. Seuils de tempratures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Frvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biolgico. Biologische Gefhrdung. Rishio biologico. Risco biolgico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instrues de utilizao. Se brugsanvisning. Se bruksanvisningen. Ls instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif mdical de diagnostic in vitro. Producto sanitario para diagnstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos mdicos para diagnstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atensin. Achtung. Attenzione. Ateno. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend fr 96 Tests. Contenuto sufficiente per 96 test. Inneholder tilstrkkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehller tillrckligt fr 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformment  la directive europenne 98/79/CE relative aux dispositifs mdicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformit alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rdets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. verensstmmer med direktiv 98/79/EG fr medicintekniske produkter.

Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigénio. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Fortynning. Diluent. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Conjugado. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x conc. Soluzione di lavaggio 30x conc. Solução de lavagem concentrada 30x. Vaskebuffer 30x koncentreret.. Vaskebuffer, konsentret 30 ganger. Tvättbuffert 30x conc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
SOLN STOP	Stop solution. Solution d'Arrêt. Solución de Parada. Stopplösung. Soluzione bloccante. Solução de paragem. Stoppløsning. Stoppløsning. Stoppløsning
CAL X	Calibrator. Etalon. Calibrador. Calibratore. Calibrador. Kalibrator. Kalibrator. Kalibrator. Kalibrator.
CONTROL X	Control. Contrôle. Controllo. Kontrolle. Controllo. Controllo. Kontroll. Kontroll. Kontroll.

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References:

1. **Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): Manual of biological markers of disease. Dordrecht, Kluwer Academic Publishers; 1993, A9:1-14.
2. **Rasmussen N, Sjolín C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). *J Immunol Methods* 1990,127:139-145.
3. **Wieslander J.** How are anti neutrophil cytoplasmic autoantibodies detected? *Am J Kidney Dis* 1991,18:154-158.
4. **Jenne DE, Tschopp J, Ludemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. *Nature* 1990, 346:520.
5. **Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
6. **Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
7. **Jennette C, Falk R.** Small vessels vasculitis. *N Engl J Med* 1997, 337:1512-1525.
8. **Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. *Kidney Int* 1995, 48: 844-850.
9. **Cohen Tervaert JW, Huitema MG, Hene RJ, Sluiter WJ, The T, van der Hem GK et al.** Prevention of relapses in Wegener's granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre. *Lancet* 1990; 336:709-711.
10. **Savige J, Davis D, Falk DJ, Jennette JC, Wiik A.** ANCA and associated diseases; A review of clinical and laboratory features. *Kidney Int* 2000, 57:846-862.
11. **Savige J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. *Am J Clin. Pathol* 1999, 111:507-513.
12. **Csernok E, Ahlquist D, Ullrich S, Gross W.L.** A critical evaluation of commercial immunoassays for antineutrophil cytoplasmic antibodies directed against proteinase 3 and myeloperoxidase in Wegener's granulomatosis and microscopic polyangiitis. *Rheumatol* 2002, 41:1313-1317.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Explicação dos símbolos Forklaring til symboler. Symboler som brukes på etiketter. Förklaringar till symboler.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partinummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Utløpsdato. Utløpsdato. Använd före.
	Temperature limit. Seuils de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbeverages ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Inneholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigénio. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Fortynning. Diluent. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Conjugado. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x conc. Soluzione di lavaggio 30x conc. Solução de lavagem concentrada 30x. Vaskebuffer 30x koncentreret.. Vaskebuffer, konsentret 30 ganger. Tvättbuffert 30x conc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
SOLN STOP	Stop solution. Solution d'Arrêt. Solución de Parada. Stopplösung. Soluzione bloccante. Solução de paragem. Stopløsning. Stoppløsning. Stopplösning
CAL X	Calibrator. Etalon. Calibrador. Calibratore. Calibrador. Kalibrator. Kalibrator. Kalibrator. Kalibrator.
CONTROL X	Control. Contrôle. Controllo. Kontrolle. Controllo. Controllo. Kontroll. Kontroll. Kontroll.

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References:

1. **Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): Manual of biological markers of disease. Dordrecht, Kluwer Academic Publishers; 1993, A9:1-14.
2. **Rasmussen N, Sjolin C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). J Immunol Methods 1990,127:139-145.
3. **Wieslander J.** How are anti neutrophil cytoplasmic autoantibodies detected? Am J Kidney Dis 1991,18:154-158.
4. **Jenne DE, Tschopp J, Ludemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. Nature 1990, 346:520.
5. **Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
6. **Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
7. **Jennette C, Falk R.** Small vessels vasculitis. N Engl J Med 1997, 337:1512-1525.
8. **Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. Kidney Int 1995, 48: 844-850.
9. **Cohen Tervaert JW, Huitema MG, Hene RJ, Sluiter WJ, The T, van der Hem GK et al.** Prevention of relapses in Wegener's granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre. Lancet 1990; 336:709-711.
10. **Savige J, Davis D, Falk DJ, Jennette JC, Wiik A.** ANCA and associated diseases; A review of clinical and laboratory features. Kidney Int 2000, 57:846-862.
11. **Savige J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. Am J Clin. Pathol 1999, 111:507-513.
12. **Csernok E, Ahlquist D, Ullrich S, Gross W.L.** A critical evaluation of commercial immunoassays for antineutrophil cytoplasmic antibodies directed against proteinase 3 and myeloperoxidase in Wegener's granulomatosis and microscopic polyangiitis. Rheumatol 2002, 41:1313-1317.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Explicação dos símbolos Forklaring til symboler. Symboler som brukes på etiketter. Förklaringar till symboler.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partinummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Utløpsdato. Utløpsdato. Använd före.
	Temperature limit. Seuils de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbeverages ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
 96	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Inneholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigénio. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Fortynning. Diluent. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Conjugado. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x conc. Soluzione di lavaggio 30x conc. Solução de lavagem concentrada 30x. Vaskebuffer 30x koncentreret.. Vaskebuffer, konsentret 30 ganger. Tvättbuffert 30x conc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
SOLN STOP	Stop solution. Solution d'Arrêt. Solución de Parada. Stopplösung. Soluzione bloccante. Solução de paragem. Stopløsning. Stoppløsning. Stopplösning
CAL X	Calibrator. Etalon. Calibrador. Calibratore. Calibrador. Kalibrator. Kalibrator. Kalibrator. Kalibrator.
CONTROL X	Control. Contrôle. Controllo. Kontrolle. Controllo. Controllo. Kontroll. Kontroll. Kontroll.

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Advarsel

Indeholder ProClin 300:
Blanding af: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

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SVAR LIFE SCIENCE AB