Instruction

**WIESLAB® PR3-ANCA**

Enzyme immunoassay for detection of autoantibodies  
against proteinase 3

Microtitration 96 wells  
Store the kit at +2-8° C  
For in vitro diagnostic use only

Document No. LABEL-DOC-0045 2.0

<b>WIESLAB® PR3-ANCA</b>	
English:	page ..... 2
Français:	page ..... 12
Español:	página ..... 19
Deutsch:	Seite ..... 26
Italiano:	pagina ..... 33
Português:	página ..... 40
Dansk:	side..... 47
Norsk:	side..... 54
Svenska:	sida..... 61

PR3 IU

**INTENDED USE**

The Wieslab® PR3-ANCA test kit is an enzyme-linked immunosorbent assay (ELISA) for detection and quantitation of IgG antibodies to proteinase 3 (PR3) in human sera. The assay is used to detect antibodies in a single serum specimen. The results of the assay are to be used as an aid to the diagnosis of Wegener's granulomatosis. The analysis should be performed by trained laboratory professionals.

FOR IN VITRO DIAGNOSTIC USE.

**Summary and explanation**

ANCAs (anti-neutrophil cytoplasmic antibodies) are a family of autoantibodies related to vasculitis and inflammatory disorders. Since 1985, when c-ANCA was shown to be related to Wegener's granulomatosis (WG) (1), interest in ANCAs has steadily increased, and today these antibodies are considered to be major diagnostic tools for the diagnosis and follow up of systemic vasculitis (13, 14). The first method to detect ANCA was indirect immuno-fluorescence (IIF) performed on ethanol fixed granulocytes (2). This method yields two patterns, a cytoplasmic staining of the granulocyte denoting the presence of c-ANCAs, and a perinuclear staining denoting the presence of p-ANCAs. IIF was followed by ELISAs using the purified proteins (3,4).

The granulocyte is full of granules each with many different proteins. It was early shown that antibodies from systemic vasculitis patients bind to the alpha fraction containing the azurophil granules. The most important proteins were PR3 (5, 6) and myeloperoxidase (MPO) (7, 8). PR3 is a serine protease with a molecular weight of 29kD, and MPO is a dimer with a molecular weight of 140kD. Thus antibodies to proteinase 3 are termed PR3-ANCA, and antibodies to myeloperoxidase are termed MPO-ANCA.

Approximately 80-90% of WG patients manifest PR3-ANCA and 5-15% MPO-ANCA. One category of vasculitis is microscopic polyangiitis (MP). Most patients with active MP are characterised by positive ANCA test results, MPO-ANCA being more frequent than PR3-ANCA (9, 10, 11, 12).

An international workgroup has developed an international standard for PR3-ANCA serology. The Wieslab® PR3-ANCA is standardized against the AF-CDC international standard, code IS2721 (Human Reference Serum 16).

**Principle of the Wieslab® PR3-ANCA**

The wells of the microtitre plate are coated with purified proteinase 3. During the first incubation, specific antibodies in diluted serum, will bind to the antigen coating.

The wells are then washed to remove unbound antibodies and other components.

A conjugate of alkaline phosphatase-labelled antibodies to human IgG binds to the antibodies in the wells in this second incubation.

After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. The amount of bound antibodies correlates to the colour intensity and is measured in terms of absorbance (optical density (OD)). The absorbance is then calculated against a calibrator curve and the results are given in IU/mL adapted to the AF-CDC international standard for PR3.

**Warnings and precautions**

- For in vitro diagnostic use.
- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommend that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The concentrations of anti-PR3 in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity.

- Material safety data sheets for all hazardous components contained in this kit are available on request from Svar Life Science.

BUF	WASH	30X
DIL		
CONJ		
CAL		

CONTROL	+
CONTROL	-
SOLN	STOP
SUBS	pNPP

### Warning

Contains ProClin 300:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

H317:	May cause an allergic skin reaction.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P302+352:	IF ON SKIN: Wash with plenty of soap and water.
P333+313:	If skin irritation or rash occurs: Get medical advice/attention.

### Specimen collection

The PR3-ANCA assay is for use with serum. Handle as if capable of transmitting infectious agents. Avoid using sera which are icteric, lipemic and hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used.

Store serum between 2-8° C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20° C or colder. Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.

The CLSI provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

### Kit components and storage of reagents

- One frame with 96 wells (red coloured) coated with proteinase 3 with lid sealed in a foil pack with a dry pack.
- 1.5 mL negative control (NC) containing human serum in diluent.
- 1.5 mL positive control (PC) containing human serum in diluent.
- 13 mL conjugate containing alkaline phosphatase-labelled antibodies to human IgG (blue colour).
- 2 x 32 mL Diluent (Dil) containing PBS (red colour).
- 13 mL Substrate pNPP.
- 13 mL Stop solution
- 30 mL wash solution 30x concentrated.
- Six calibrators (five calibrators, Cal 2-6, containing human serum) in diluent. 1.5 mL Cal 1 = 0 IU/mL, 1.5 mL Cal 2 = 2 IU/mL, 1.5 mL Cal 3 = 10 IU/mL, 1.5 mL Cal 4 = 30 IU/mL, 1.5 mL Cal 5 = 100 IU/mL, 1.5 mL Cal 6 = 200 IU/mL.

All reagents in the kit are ready for use except wash solution and should be stored at 2-8° C. Remove only the number of strips needed for testing, resealing the aluminium package carefully.

### Materials or equipment required but not provided

- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- Washer for strips, absorbent tissue, tubes and a timer.

**PROCEDURE**

All solutions should be used at room temperature. Incubate all steps at room temperature (18-25° C) with lid. Incubate serum for 60 minutes, conjugate for 30 minutes and substrate for 30 minutes.

**Preparation of washing solution**

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37° C water bath until the crystals have dissolved before dilution of wash solution. Dilute 10 mL of the 30x concentrated wash solution in 290 mL distilled water. When stored at 2-8° C, the diluted wash solution is stable until the date of expiration of the kit.

**Dilution of serum and incubation**

Dilute the patient's serum 1/100 with diluent (990 µL diluent + 10 µL serum).

Pipet 100µL/well in duplicate of Calibrator 1, 2, 3, 4, 5, 6, PC, NC and diluted patient's serum (P) according to the diagram below. Incubate for 60 minutes.

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	Cal 1	Cal 5	P1									
<b>B</b>	Cal 1	Cal 5	P1									
<b>C</b>	Cal 2	Cal 6	P2									
<b>D</b>	Cal 2	Cal 6	P2									
<b>E</b>	Cal 3	PC	etc									
<b>F</b>	Cal 3	PC										
<b>G</b>	Cal 4	NC										
<b>H</b>	Cal 4	NC										

**After serum incubation**

Wash 3 times with 300 µL washing solution / well, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

**Adding conjugate**

Add 100 µL conjugate to each well. Incubate for 30 minutes.

**After conjugate incubation**

Wash as before.

**Adding substrate solution**

Add 100 µL substrate pNPP to each well, incubate for 30 minutes.

**Adding stop solution**


Add 100 µL stop solution to each well and read the absorbance at 405 nm on a microtiterplate reader within 2 hours.

**Calculations**

Construct a calibrator curve by plotting the OD against the IU/mL values of the 6 calibrators. The six calibrators provided have been titrated to match the AF-CDC standard resulting in the values of 0 IU/mL for calibrator 1, 2 IU/mL for calibrator 2, 10 IU/mL for calibrator 3, 30 IU/mL for calibrator 4, 100 IU/mL for calibrator 5 and 200 IU/mL for calibrator 6. Read the IU/mL value of the patient from the constructed curve. Values greater than 200 should be reported as >200, or reassay them with a higher dilution. The Wieslab® PR3-ANCA is standardized against the AF-CDC international standard.

<u>Example:</u>	<u>Calibrator</u>	<u>IU/mL</u>	<u>Absorbance</u>
	1	0	0.102
	2	2	0.192
	3	10	0.509
	4	30	0.902
	5	100	1.640
	6	200	1.906

A sample with an absorbance value of 0.886 will read on the X-axis as having 27 IU/mL of PR3-ANCA. In this example a 4 parameter logistic curve fit has been applied.

**Important:** The curve is an example and should not be used for actual patient interpretation.

### Quality Control

The OD for calibrator 1 should be < 0.2.

The OD for calibrator 6 should be > 1.0

The value for the positive and negative controls see lot certificate.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cut-off. It is recommended to assay an additional control at the assay cut-off. As the positive control is ready-to-use it does not indicate an eventual dilution error by the user. It is recommended to use an internal control for this purpose.

If any of the values are not within their respective ranges, the test should be considered invalid and the test should be repeated. Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organisations. Refer to CLSI C24-A for guidance on appropriate QC practices.

### Interpretation of results

< 4 IU/mL = **Negative**

4-6 IU/mL = **Equivocal**; Retest, if still equivocal retest by an alternative method or test a new sample

> 6 IU/mL = **Positive**

### Limitations

The individual patient's antibody titre can not be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardisation of results.

The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used as an adjunct to clinical symptoms and the results of other available tests.

Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive for PR3 antibodies with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies.

Immunosuppressive therapy should not be started on basis of a positive ANCA result. Initiation or changes in treatment should not be based on changes in ANCA concentration alone, but rather on careful clinical observation.

### Expected results

PR3-ANCA are rarely found in normal healthy individuals. The PR3-ANCA was tested with 120 normal sera, of which all were found to be negative. One new patient with PR3-ANCA is expected per 100000 individuals per year. Around 10% of patients with WG are negative in both IIF and ELISA. The PR3-ANCA was tested with 83 WG sera, 64 were found to be positive. PR3-ANCA are found in approximately one half of the sera of MP (microscopic polyangiitis) patients. The PR3-ANCA was tested with 55 MP sera, 21 were found to be positive.

## Performance characteristics

**Table 1. Clinical sensitivity and specificity.** A total of 353 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total number	Negative <4 IU/mL	Equivocal 4-6 IU/mL	Positive >6 IU/mL
Blood donors	120	120	0	0
WG	83	12	6	65
MP	55	30	4	21
SLE	14	13	1	0
Others	26	25	0	1
Anti-GBM	55	52	2	1

WG = Wegener's granulomatosis,

MP = microscopic polyangiitis

SLE = systemic lupus erythematosus

Others = RA, Sjögren's syndrome

GBM = glomerular basement membrane

### Clinical sensitivity (equivocal samples excluded from calculations)

WG =  $65/77 = 84.4\%$       95% CI = 74.4-91.7 %

MP =  $21/51 = 41.2\%$       95% CI = 27.6-55.8 %

### Clinical specificity (equivocal samples excluded from calculations)

GBM =  $52/53 = 98.1\%$       95% CI = 89.9-100 %

SLE =  $13/13 = 100\%$       95% CI = 75.3-100 %

Others =  $25/26 = 96.2\%$       95% CI = 80.4-99.9 %

Donors =  $120/120 = 100\%$       95% CI = 97.0-100 %

The 95% confidence interval (CI) was calculated using the exact method.

**Table 2. Agreement between the Wieslab® PR3-ANCA and an alternative ELISA.** A total of 353 frozen retrospective sera were assayed. The following table summarises the results.

Wieslab® PR3-ANCA					
Alternative ELISA		Positive	Equivocal	Negative	Total
	Positive	77	0	0	77
	Equivocal	10	4	2	16
	Negative	1	9	250	260
	Total	88	13	252	353

Sera falling in the equivocal range were excluded from the following calculations.

Positive Percent Agreement =  $77/77 = 100\%$       95% CI = 95.3-100 %

Negative Percent Agreement =  $250/251 = 99.6\%$       95% CI = 97.8-100 %

Overall Percent Agreement =  $327/328 = 99.7\%$       95% CI = 98.3-100 %

The 95% confidence interval (CI) was calculated using the exact method.

**Table 3. Batch to batch variation** was determined by testing six different samples in eight replicates on three different batches.

	1	2	3	4	5	6
Mean value (IU/mL)	8	12	17	20	23	65
SD	0.39	0.02	0.27	0.99	1.86	10.70
% CV	5	0.1	2	5	8	16

**Table 4. Inter-assay precision** was determined by testing six different samples in eight replicates at three separate occasions.

	1	2	3	4	5	6
Mean value (IU/mL)	8	13	15	20	35	60
SD	0.41	0.75	3.16	1.47	6.12	11.69
% CV	5	6	21	7	17	19

**Table 5. Intra-assay precision** was determined by testing six different samples in eight replicates at one occasion.

	1	2	3	4	5	6
Mean value (IU/mL)	7	12	16	19	21	58
SD	0.45	0.75	1.74	0.92	0.75	6.62
% CV	6	6	11	5	4	11

**Table 6. Dilution recovery** was determined by testing five serial dilutions for three different samples.

Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution Corrected % Recovery
1	1/100	62	62	100
	1/200	29	31	94
	1/400	14	16	90
	1/800	7	8	90
	1/1600	3	4	77
Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution corrected % Recovery
2	1/100	118	118	100
	1/200	54	59	92
	1/400	26	30	88
	1/800	12	15	81
	1/1600	6	7	81
Sample	Dilution	Mean Measured Concentration (IU/mL)	Calculated Concentration (IU/mL)	Dilution corrected % Recovery
3	1/100	139	139	100
	1/200	58	70	83
	1/400	27	35	78
	1/800	14	17	81
	1/1600	7	9	81

Since patient sera contain heterogeneous antibody populations some samples may exhibit non-linearity, especially at very high sample dilutions.

**Troubleshooting**

<b>Problem:</b>	<b>Possible causes:</b>	<b>Solution:</b>
Control values out of range.	<ol style="list-style-type: none"> <li>1. Incorrect temperature, timing or pipetting; reagents not mixed.</li> <li>2. Cross contamination of controls.</li> <li>3. Improper dilution.</li> <li>4. Optical pathway not clean.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test.</li> <li>2. Pipette carefully.</li> <li>3. Repeat test.</li> <li>4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.</li> </ol>
All test results negative.	<ol style="list-style-type: none"> <li>1. One or more reagents not added, or added in wrong sequence.</li> <li>2. Antigen coated plate inactive.</li> </ol>	<ol style="list-style-type: none"> <li>1. Recheck procedure. Check for unused reagent. Repeat test.</li> <li>2. Check for obvious moisture in unused wells. Wipe bottom and reread.</li> </ol>
All test results yellow.	<ol style="list-style-type: none"> <li>1. Contaminated buffers or reagents.</li> <li>2. Washing solution contaminated.</li> <li>3. Improper dilution of serum.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check all solutions for turbidity.</li> <li>2. Use clean container. Check quality of water solution used to prepare.</li> <li>3. Repeat test.</li> </ol>
Poor precision.	<ol style="list-style-type: none"> <li>1. Pipette delivery CV greater than 5%.</li> <li>2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature.</li> <li>3. Reagent addition taking too long; inconsistency in timing intervals.</li> <li>4. Optical pathway not clean.</li> <li>5. Washing not consistent; trapped bubbles; washing solution left in the wells.</li> <li>6. Improper pipetting.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check calibration of pipette. Use reproducible technique.</li> <li>2. Mix all reagents gently but thoroughly and equilibrate to room temperature.</li> <li>3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time.</li> <li>4. Check for air bubbles in the wells. Wipe bottom and reread.</li> <li>5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue.</li> <li>6. Avoid air bubbles in pipette tips.</li> </ol>

**References:**

1. **van der Woude FJ, et al.** Autoantibodies against neutrophils and monocytes: Tool for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet* 1985, i 425-429.
2. **Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): *Manual of biological markers of disease.* Kluwer Academic Publishers 1993, A9, 1-14.
3. **Rasmussen N, Sjölin C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). *J Immunol Methods* 1990, 127, 139-145.
4. **Segelmark M, et al.** How and why should we detect ANCA? *Clin Exp Rheumatol* 2000, 18, 629-635.
5. **Jenne DE, Tschopp J, Lüdemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. *Nature* 1990, 346, 520.
6. **Wieslander J, Wik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease.* Kluwer Academic Publishers 1994, 1-9.
7. **Falk RJ, Jennette JC.** Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. *N Engl J Med* 1988, 318, 1651-1657.
8. **Wieslander J, Wik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease.* Kluwer Academic Publishers 1994, 1-9.
9. **Kamech L, Harper L, Savage C.** ANCA - Positive Vasculitis. I. *Am Soc Nephrol* 2002 43 1953 - 1960.
10. **Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. *Kidney Int* 1995, 48, 844-850.
11. **Boomsma MM, Stegeman CA, van der Leij MJ, et al.** Prediction of relapses in Wegener's granulomatosis by measurement of ANCA levels: a prospective study. *Arthritis Rheum*, 2000, 43 2025 - 2033.
12. **Savige J, Davies D, Falk RJ, Jennette JC, and Wik A.** ANCA and associated diseases: A review of the clinical and laboratory features. *Kidney Int* 2000, 57, 846-862.
13. **Savige J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. *Am J Clin Pathol* 1999, 111, 507-513.
14. **Won H, et al.** Serial ANCA titers. Useful tool for the prevention of relapses in ANCA associated vasculitis. *Kidney Int*, 2003, 63, 1079-1085.
15. **Neumann, I, et al.** Glomerular immune deposits are associated with increased proteinuria in patients with ANCA-associated crescentic nephritis. *Nephrol Dial Transplant* 2003, 18, 524-531.
16. **Csernok E, Ahlquist D, Ullrich S, Gross W.L.** A critical evaluation of commercial immunoassays for antineutrophil cytoplasmic antibodies directed against proteinase 3 and myeloperoxidase in Wegener's granulomatosis and microscopic polyangiitis. *Rheumatology* 2002;41:1313-1317.

**Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Explicação dos símbolos Forklaring til symboler. Symboler som brukes på etiketter. Förklaringar till symboler.**

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partinummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Utløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuil de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

<b>Ag</b>	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigénio. Antigen. Antigen. Antigen.
<b>DIL</b>	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Fortynning. Diluent. Spädningsbuffert.
<b>CONJ</b>	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Conjugado. Konjugat. Konjugat. Konjugat.
<b>BUF</b> <b>WASH</b> <b>30X</b>	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Solução de lavagem concentrada 30x. Vaskebuffer 30x koncentreret.. Vaskebuffer, konsentrert 30 ganger. Tvättbuffert 30x konc.
<b>SUBS</b> <b>pNPP</b>	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
<b>SOLN</b> <b>STOP</b>	Stop solution. Solution d'Arrêt. Solución de Parada. Stopplösung. Soluzione bloccante. Solução de paragem. Stoppløsning. Stoppløsning. Stoppløsning.
<b>CAL</b> <b>X</b>	Calibrator. Etalon. Calibrador. Calibratore. Calibrador. Kalibrator. Kalibrator. Kalibrator. Kalibrator.
<b>CONTROL</b> <b>X</b>	Control. Contrôle. Controllo. Kontrolle. Controllo. Controllo. Kontroll. Kontroll. Kontroll.

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**References:**

1. **van der Woude FJ, et al.** Autoantibodies against neutrophils and monocytes: Tool for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet* 1985, i 425-429.
2. **Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): *Manual of biological markers of disease*. Kluwer Academic Publishers 1993, A9, 1-14.
3. **Rasmussen N, Sjölin C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). *J Immunol Methods* 1990, 127, 139-145.
4. **Segelmark M, et al.** How and why should we detect ANCA? *Clin Exp Rheumatol* 2000, 18, 629-635.
5. **Jenne DE, Tschopp J, Lüdemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. *Nature* 1990, 346, 520.
6. **Wieslander J, Wik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Kluwer Academic Publishers 1994, 1-9.
7. **Falk RJ, Jennette JC.** Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. *N Engl J Med* 1988, 318, 1651-1657.
8. **Wieslander J, Wik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Kluwer Academic Publishers 1994, 1-9.
9. **Kamech L, Harper L, Savage C.** ANCA - Positive Vasculitis. I. *Am Soc Nephrol* 2002 43 1953 - 1960.
10. **Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. *Kidney Int* 1995, 48, 844-850.
11. **Boomsma MM, Stegeman CA, van der Leij MJ, et al.** Prediction of relapses in Wegener's granulomatosis by measurement of ANCA levels: a prospective study. *Arthritis Rheum*, 2000, 43 2025 - 2033.
12. **Savige J, Davies D, Falk RJ, Jennette JC, and Wik A.** ANCA and associated diseases: A review of the clinical and laboratory features. *Kidney Int* 2000, 57, 846-862.
13. **Savige J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. *Am J Clin Pathol* 1999, 111, 507-513.
14. **Won H, et al.** Serial ANCA titers. Useful tool for the prevention of relapses in ANCA associated vasculitis. *Kidney Int*, 2003, 63, 1079-1085.
15. **Neumann, I, et al.** Glomerular immune deposits are associated with increased proteinuria in patients with ANCA-associated crescentic nephritis. *Nephrol Dial Transplant* 2003, 18, 524-531.
16. **Csernok E, Ahlquist D, Ullrich S, Gross W.L.** A critical evaluation of commercial immunoassays for antineutrophil cytoplasmic antibodies directed against proteinase 3 and myeloperoxidase in Wegener's granulomatosis and microscopic polyangiitis. *Rheumatology* 2002;41:1313-1317.

**Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Explicação dos símbolos Forklaring til symboler. Symboler som brukes på etiketter. Förklaringar till symboler.**

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partinummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Utløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuil de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

<b>Ag</b>	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigénio. Antigen. Antigen. Antigen.
<b>DIL</b>	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Fortynning. Diluent. Spädningsbuffert.
<b>CONJ</b>	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Conjugado. Konjugat. Konjugat. Konjugat.
<b>BUF</b> <b>WASH</b> <b>30X</b>	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Solução de lavagem concentrada 30x. Vaskebuffer 30x koncentreret.. Vaskebuffer, konsentrert 30 ganger. Tvättbuffert 30x konc.
<b>SUBS</b> <b>pNPP</b>	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
<b>SOLN</b> <b>STOP</b>	Stop solution. Solution d'Arrêt. Solución de Parada. Stopplösung. Soluzione bloccante. Solução de paragem. Stoppløsning. Stoppløsning. Stoppløsning.
<b>CAL</b> <b>X</b>	Calibrator. Etalon. Calibrador. Calibratore. Calibrador. Kalibrator. Kalibrator. Kalibrator. Kalibrator.
<b>CONTROL</b> <b>X</b>	Control. Contrôle. Controllo. Kontrolle. Controllo. Controllo. Kontroll. Kontroll. Kontroll.

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**SVAR LIFE SCIENCE AB**